



Labpak


Everything For
Your Laboratory

Presspak

SPECIAL OFFERS

SPRING 2009

New from Mettler!



New FiveGo™ pH - FG2 Food Kit

The FG2 portable pH meter provides rapid and reliable results wherever you are. This new portable kit version contains the easy and compact pH meter FiveGo™ in combination with the new LE427 pH puncture electrode which is tailor made for portable pH measurements in the food industry and for related puncture applications.

The new LE427 puncture pH electrode with its spear shaped tip makes it easy to penetrate solid or semi-solid samples such as food and soil and comes with a 1.2m cable

- IP54 Protection against dust & water
- Fast changeover between pH and mV
- Data storage for 30 measurements
- Automatic calibration



FG2-Foodkit

**Normally £220.00
NOW ONLY £175.00**

Order using item ref 51302906. Comes complete with Instrument, electrode clip, wrist strap, manual, LE427 electrode, IP54 connection caps and buffer sachets.

Food Measurement

Mettler-Toledo's successful and competitively priced refractometers, the Quick-Brix™ 60 and Quick-Brix™ 90 simplify and speed up processes in the food industry. Both models are designed for ease of use with:-

- Simple operation
- Compact size - fits in every pocket
- Fast results
- Easy to clean design

Ingredients testing saves money and improves food product quality. Quick-Brix™ is easy to use, accurate and reliable. Wipe the prism with a tissue for cleaning and calibrate with a single key stroke and water as a standard. As easy as 1,2,3.

Quick-Brix™ 60

- 0-60% Brix
- Resolution of 0.1% Brix, Accuracy of 0.2%.

Quick-Brix™ 90

- 0-90% Brix
- Resolution of 0.1% Brix, Accuracy of 0.1%



The Quick-Brix™

**Free Thermal Mug with
any item off this page!**

QUICKBRIX 60 ONLY £205. QUICKBRIX 90 ONLY £435

**Unit 6
Ketterer Court
Jackson Street
St Helens
Merseyside
WA9 3AH**

Offers valid until end of June 2009

For more information on any of these offers or to arrange a quotation, please do not hesitate to contact us by phone, fax or e-mail.

**Phone: 0870 034 2055
Fax: 0870 034 2056
E-mail: sales@prlabs.co.uk
Web: www.prlabs.co.uk**