



Labpak
Everything For
Your Laboratory

Presspak

SPECIAL OFFERS

WINTER 2009

BATH ARMOR™ FULL METAL

Still using water? Still using solid heat blocks?

With Bath Armor™ you don't need to.

Replace water in your water bath with these thermally efficient pellets that are easy to keep clean. It enables you to eliminate nearly all water bath cleaning and protects your water bath from drying out. It saves the need for racks, floaties and bottle weights. It means more free time, one less lab duty and you never have to wash or refill a water bath again.



New Bath Armor™

- Moisture and Gas impermeable
- High thermal conductivity
- Smooth, rounded and fluid like
- Reduces water consumption
- No more filling
- No waterborne contamination
- Washes & sterilizes easily
- Good thermal heat distribution



Replaces pre-drilled blocks in your block heaters with thermally efficient pellets that are easy to keep clean. Reduces investment in blocks and heating well can be used for heating any tube size or shape at the same time.

- One pile suits all tube sizes
- Replaces pre-drilled blocks
- More convenient to use

Ref BA-5 - £235
Ref BA-0.5 - £50

Please note—Ref BA-5—One pack required for 2, 4 & 5L baths giving fill depth of 80mm.
For an 8L bath, 1 pack gives depth of 40mm or 2 packs for twice the fill depth. BA-0.5 is for block heaters.

New for 2009 - Clifton's Thermostatic Hotplates which have been completely redesigned. Three models - small HPI-1, medium HPI-2 and large HPI-3. The hob is solid aluminium 12mm thick for excellent thermal properties and uniform heat distribution over its working area. All painted surfaces and the control panel are antibacterial & the controls have been updated with the addition of a 'Hot' warning indicator.

Accessory Sand bath kits are available as an alternative to mineral-oil baths used at higher temperatures. There's no open container of liquid to worry about and sand won't evaporate or impose temperature limits. This is especially useful when something must be heated above 250 °C (the approximate limit for an oil bath).

Digital versions are also available. Just ask for details.



The NEW HPI-2D

HPI-1 - £425
HPI-2 - £485
HPI-3 - £645

Unit 6
Ketterer Court
Jackson Street
St Helens
Merseyside
WA9 3AH

[Offers valid until end of March 2009](#)

For more information on any of these offers or to arrange quotation, please do not hesitate to contact us by phone, fax or e-mail.

Phone: 0870 034 2055
Fax: 0870 034 2056
E-mail: sales@prilabs.co.uk
Web: www.prilabs.co.uk