

## Limited Stock Offers

Available only whilst stocks lasts!



### Stuart SM4 mini magnetic stirrer

- Compact design (90x160x90mm (wxdxh))
- Speed range: 80 to 1300rpm
- Suitable for 16mm retort rod

### IKA BIG SQUID stirrers

- Speed range: 0 to 1500rpm
- Stirring capacity: 1500ml (H<sub>2</sub>O).
- Plate diameter: 160mm
- Dimensions: 180x180x48mm (wxdxh)



### Ika lab disc Magnetic Stirrers

Ultra-flat, compact magnetic stirrer featuring modern magnet coil technology and a wear-free drive (no moving parts). The lab disc can reverse direction of rotation automatically every 30 seconds to ensure better mixing.

### VWR Collection Shaker Digital Lab Dancer



Economic, compact test tube shaker with touch function with 6 different speed settings at the touch of a button:  
Level 1 approx. 1000 rpm to level 6 approx. 2800 rpm The selected speed level is shown on an LED display.

- The lab dancer can be used with all small vessels of up to 30 mm in diameter, e.g. centrifuge tubes, Eppendorf vessels and test tubes up to 50 ml
- The upper casing and the test tube surface are made from inert plastic

Description	Pk.	Cat. No.	List Price £	Promo Price £
STIRRER MAGNETIC MINI SM4	1 * 1 items	442-0308	137.00	80.00
IKA BIG SQUID stirrers - STAR	1 * 1 items	442-0034	135.00	80.00
IKA BIG SQUID stirrers - FROGGY	1 * 1 items	442-0035	135.00	80.00
IKA BIG SQUID stirrers - OCEAN	1 * 1 items	442-0037	135.00	80.00
IKA BIG SQUID stirrers - IKAMALEON	1 * 1 items	442-0038	135.00	80.00
Lab disc package tennis player (uk plug) with exclusive rosenthal city-cup and t-shirt: new ultra-flat magnetic stirrer	1 * 1 items	442-0209	120.00	70.00
Lab disc magnetic stirrer air-mail design	1 * 1 items	442-0289	98.50	60.00
Lab disc magnetic stirrer mail-box design	1 * 1 items	442-0288	98.50	60.00
Lab disc magnetic stirrer stone washed design	1 * 1 ITEMS	442-0290	98.50	60.00
Shaker test tube lab dancer digital	1 * 1 items	444-2338	137.00	80.00

**Order from:**  
**P&R Labpak Ltd.**  
 Unit 6 Ketterer Court  
 Jackson Street  
 St Helens  
 Merseyside  
 WA9 3AH

T: 0870 034 2055  
 F: 0870 034 2056  
 E: sales@prlabs.co.uk  
 W: www.prlabs.co.uk