



**Labpak**

Everything For  
Your Laboratory

SPECIAL OFFERS

SPRING 2009

# Presspak


## New from Mettler!

### New SevenGo Duo Pro™



The NEW SG68

**£20 Argos Voucher  
with each unit bought**



The new professional multi-parameter instruments cover the important parameters pH, ion concentration, conductivity and dissolved oxygen. Handy features like the automatic air pressure compensation, the helpful ion mode and the wireless infrared communication port save time and guarantee correct results.

Mistakes are easily avoided with the innovative Intelligent Sensor Management ISM®. Correct identification of the pH, conductivity or dissolved oxygen sensor is guaranteed upon connection to the meter. The sensor type and serial number is automatically recognised.

The calibration data including date and time is stored on the sensor. This eliminates mistakes when transferring the sensor to another instrument.

On offer is the SG68-FK2 SevenGo Duo Pro™ pH/ion/dissolved oxygen (2m) field kit (ref 51302612). Complete with IP67 electrodes, field case, electrode clips, buffer sachets and rubber holster.

**Normally £1436. NOW ONLY £1149**

### New SevenGo Duo™



The NEW SG23

The new SevenGo Duo™ are a mid-range of dual channel portable meters for pH and conductivity measurements. They are ergonomic in design and have a simple user interface for efficiency in the field and come with Intelligent Sensor Management® for extra security. Other features include:-

- Data storage for 99 GLP data records
- IP67 ingress protection
- Automatic end point recognition

**£20 Argos Voucher  
with each unit bought**

On offer is the SG23-FK2 SevenGo Duo™ pH/conductivity (2m) field kit (ref 51302602). Complete with IP67 electrodes, field case, electrode clips, buffer sachets and rubber holster and more.

**Normally £892. NOW ONLY £715**

**Unit 6  
Ketterer Court  
Jackson Street  
St Helens  
Merseyside  
WA9 3AH**

Offers valid until end of June 2009

For more information on any of these offers or to arrange a quotation, please do not hesitate to contact us by phone, fax or e-mail.

**Phone: 0870 034 2055  
Fax: 0870 034 2056  
E-mail: sales@prilabs.co.uk  
Web: www.prilabs.co.uk**