

Sealpette Pro, Variable Volume, Single Channel pipettor

Sealpette Pro, Variable Volume, Single Pipettor



Fully autoclavable
 Very low pipetting forces needed
 Light weight
 Tip cone filters and tip ejection by pressing the operating button
 Thermal insulation of internal components
 Small increments
 Fast and light volume adjustment
 For both right- and left-handed use
 Clear and easy-to-read display
 Personal identification with an ID tag
 Convenient and stylish pipettor shelf hanger
 Easy maintenance with one multi-purpose tool:

- calibration (standard hexagon key)
 - opening tool for the ID window
 - opening tool for tubes
- 3 year warranty

Accurate and Precise
 To ensure accuracy and precision these pipettors are individually calibrated and supplied with a quality control certificate.

Fully Autoclavable
 The Sealpette Pro pipettor is fully autoclavable (20 minutes at 121°C) without the need for re-calibration.
 These pipettors meet appropriate DIN/ISO standards. Each instrument passes stringent control procedures in accordance with ISO/DIS 8655 and bears its own QC certificate.

**Great offer
 £100!**

Description	Cat. No.	Increment, µl	Test Volume, µl	Inaccuracy (±%)	Imprecision (±%)	List Price £	Promo Price £
Sealpette Pro 0.1-3ul variable	480-044	0.002	3/1.5/0.3	1.30/2.40/10.00	0.8/1.6/6.0	171.00	100.00
Sealpette Pro 0.5-10ul variable	480-095	0.01	10/5/1/0.5	1.00/1.50/2.50/5.00	0.60/1.00/1.50/4.00	171.00	100.00
Sealpette Pro 2-20ul variable	480-039	0.02	20/10/2	0.90/1.20/3.00	0.40/1.00/2.00	161.00	100.00
Sealpette Pro 10-100ul variable	480-096	0.1	100/50/10	0.80/1.00/2.00	0.15/0.40/1.00	161.00	100.00
Sealpette Pro 20-200ul variable	480-097	0.2	200/100/20	0.60/0.80/2.00	0.15/0.30/0.80	161.00	100.00
Sealpette Pro 100-1000ul variable	480-098	1	1000/500/100	0.60/0.70/1.50	0.20/0.25/0.70	161.00	100.00
Sealpette Pro 0.1-3ul variable	480-099	10	5000/2500/500	0.50/0.60/1.50	0.15/0.30/0.60	161.00	100.00

Order from:
P&R Labpak Ltd.
 Unit 6 Ketterer Court
 Jackson Street
 St Helens
 Merseyside
 WA9 3AH

T: 0870 034 2055
 F: 0870 034 2056
 E: sales@prlabs.co.uk