



Labpak

Everything For
Your Laboratory

Presspak

SPECIAL OFFERS

AUTUMN 2009

JENWAY

**20%
OFF!**

Save an amazing 20% on the two popular Jenway bench pH Meters below and the 6305 UV/Visible Spectrophotometer!

The 3510 is a versatile, easy to use pH, mV and temperature meter ideal for routine analysis.

With 3dp resolution and a choice of 3 calibration points the meter provides the user with added flexibility. Supplied complete with glass combination electrode, electrode holder, ATC probe & buffer sachets.

- Simultaneous readout pH & temp
- 1,2 or 3 point calibration
- Auto or manual buffer recognition
- Storage of up to 32 results
- RS232 connection



The 3520 is ideal for quality control and GLP applications. Research grade specifications mean a comprehensive range of features and functions.

Its built in data logger can store, print or output readings based on a wide range of trigger options making this model suitable for monitoring and controlling pH levels.

- Ideal for quality control & GLP applications
- Storage of up to 500 results
- Multiple language options
- RS232 or infra red connection
- GLP compliant

Ref 351 001—Was £430.00 NOW £344.00

Ref 352 001—Was £608.00 NOW £486.40

Jenway 6305 UV/Visible Spectrophotometer

The 6305 is ideal for applications requiring measurements in the UV and visible wavelength ranges. Using a single, pulsed xenon lamp the 6305 offers an extended lamp life and high energy throughput in the UV and visible regions.

It is a high quality, low cost unit for spectrophotometric applications in schools and colleges, as well as in QC and testing procedures for a wide variety of services and industries.

- Wavelength range of 198 to 1000nm
- Absorbance, %T and concentration modes
- Wide range of sampling accessories
- Xenon Lamp



The JENWAY 6305

**REF 635-001
WAS £2249
NOW £1799.20**

**Unit 6
Ketterer Court
Jackson Street
St Helens
Merseyside
WA9 3AH**

Offers valid until end of Dec 2009

For more information on any of these offers or to arrange quotation, please do not hesitate to contact us by phone, fax or e-mail.

**Phone: 0870 034 2055
Fax: 0870 034 2056
E-mail: sales@prlabs.co.uk
Web: www.prlabs.co.uk**