

## VWR Collection Liquid Handling - Great Summer Offers!



612-4546

### VWR Collection Smoothie pipette filler

The VWR Smoothie pipette filler is a simple manual pipetting aid for convenient single handed use of glass or plastic pipettes from 1 up to 25 ml. The unique range selector allows to select from two dispense resolutions. One range setting delivers fast filling with optimal accuracy when used with 10 to 25 ml pipettes, the other setting is for use with pipettes from 1 to 10 ml where accuracy can sometimes be difficult to control.

- Comfortable and simple to use
- Smooth and controllable thumb action
- Silicone adapter for leakproof positioning of each pipette



612-4548

### VWR Collection Manupette pipette filler

Great for serial pipetting

The VWR Manupette pipette filler can be used with glass or plastic pipettes from 1 to 100 ml.

- Exchangeable hydrophobic membrane filter to protect the instrument from liquid penetration
- Separate buttons for aspirating / dispensing and discharging
- Silicone adapter for leakproof positioning of each pipette



612-4553

### VWR Collection Accurpette pipette controller

Fast and accurate motor-driven pipette controller.

The mode selection switch enables to choose between high or low aspirate or dispense speed. An additional gravity dispense mode is designed for use with "To Deliver" (TD) pipettes.

- Fills a 25 ml pipette in under 3 seconds on its fastest setting
- Exchangeable hydrophobic membrane filter for protection of the unit against liquid and for protection of samples against contamination
- CE marked
- Delivered with charger, 2 spare hydrophobic filters (1x 0.45 µm, 1 x 0.2 µm), bench stand/wall bracket and instruction manual



613-5287

### VWR Collection Titras Pro digital burette

Accurate high precision at a great price!

The Titras Pro digital burette combines modern design, convenient use and accuracy in one instrument. The robust outer casing does not move during the priming or dispensing sequences whilst its PTFE plunger and borosilicate glass barrel allows dispensing of aqueous solutions and some aggressive liquids.

35% off

32% off

40% off

45% off

Description	Pk.	Cat. No.	List Price £	Promo Price £
Pipette Filler Smoothie 1-25ml	1 unit	612-4546	18.40	12.00
Pipette Filler Manupette	1 unit	612-4548	36.60	25.00
Pipette Controller Accurpette, UK Plug	1 unit	612-4553	202.00	120.00
Digital Burette Titras Pro, 50ml	1 unit	613-5287	362.00	200.00

**Order from:**  
**P&R Labpak Ltd.**  
 P&R Labpak Limited  
 Unit 6 Ketterer Court  
 Jackson Street  
 St Helens  
 Merseyside  
 WA9 3AH

T: 0870 034 2055  
 F: 0870 034 2056  
 E: sales@prlabs.co.uk  
 W: www.prlabs.co.uk