



Atmosphere Generation System

Until 30
September, 2009
we are offering
great discounts
on our anaerobic
jars.



Buy 2 Get 1 FREE or Buy 1 Get 1 half price

AnaeroJar - Order code AG0025A

With a 2.5 litre capacity and supplied complete, this jar holds up to 12 plates and is compatible with the Oxoid AnaeroGen, CampyGen, CO₂Gen atmosphere generation systems, which do not require a catalyst or the addition of water.

N.B.: not suitable for use with traditional atmosphere generate systems

- 3 jars (buy 2, get 1 free) for just £270
- 2 jars (buy 1, get one half price) for £202

Anaerobic Jar - Order code HP0011A

With a 3.5 litre capacity, this jar holds up to 15 plates and is ideal for use with traditional atmosphere generation systems that do require the addition of water.

N.B.: for connection to a vacuum pump and gas cylinder, Schrader valve chucks and clips are also required

- 3 jars (buy 2, get 1 free) for just £634
- 2 jars (buy 1, get one half price) for £475

To order, ring 0870 034 2055 and quote AGS Offer.
But act now! This offer is only available until 30 September 2009.

